Current Issues in Disinfection and Sterilization

William A. Rutala, Ph.D., M.P.H. University of North Carolina (UNC) Health Care System and UNC at Chapel Hill, NC

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- Current issues
 - Critical-cleaning with washer disinfectors, Class 6 chemical indicator
 - Semicritical items (e.g., endoscopes and C. difficile spores, new AERs)
 - Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores
 - ◆ Microfibe
 - ◆ Computers-sustained antimicrobial activity, touchscreen cleaning
 - ◆ Gennicides-MRSA inactivation by disinfectants, technique
 - ◆ Green products

Disinfection and Sterilization

- EH Spaulding believed that how an object will be disinfected depended on the object's intended use.
- CRITICAL objects which enter normally sterile tissue or the vascular system or through which blood flows should be **sterile**.
- SEMICRITICAL objects that touch mucous membranes or skin that is not intact require a disinfection process (high-level disinfection [HLD]) that kills all microorganisms but high numbers of bacterial spores.
- NONCRITICAL -objects that touch only intact skin require low-level disinfection.

Cleaning

- Mechanical cleaning machines-automated equipment may increase productivity, improve cleaning effectiveness, and decrease worker exposure
 - Utensil washer-sanitizer
 - Ultrasonic cleaner
 - Washer sterilizer
 - Dishwasher
 - Washer disinfector
- Manual

Bioburden on Surgical Devices

- . Bioburden on instruments used in surgery (Nystrom, 1981)
 - 62% contaminated with <101
 - 82% contaminated with <10²
 - 91% contaminated with <103
- · Bioburden on surgical instruments (Rutala, 1997)
 - 72% contained <101
 - 86% contained <10²

Washer/Disinfector

Rutala WA, Gergen MF, Weber DJ, Unpublished results, 2007

- Five Chambers
 - Pre-wash: water/enzymatic is circulated over the load for 1 min
 - Wash: detergent wash solution (150°F) is sprayed over load for 4 min
 - Uftrasonic cleaning: basket is lowered into ultrasonic cleaning tank with detergent for 4 min
 - Thermal and lubricant rinse: hot water (180°F) is sprayed over load for 1 min; instrument milk lubricant is added to the water and is sprayed over the load
 - Drying: blower starts for 4 min and temperature in drying chamber 180F

Washer/Disinfector Removal/Inactivation of Inoculum (Exposed) on Instruments

WD Conditions	Organism	Inoculum	Log Reduction	Positives
Routine	MRSA	2.6x10 ⁷	Complete	0/8
Routine	VRE	2.6x10 ⁷	Complete	0/8
Routine	P aeruginosa	2.1x10 ⁷	Complete	0/8
Routine	M terrae	1.4x10 ⁸	7.8	2/8
Routine	GS spores	5.3x10 ⁶	4.8	11/14
No Enz/Det	VRE	2.5x10 ⁷	Complete	0/10
No Enz/Det	GS spores	8.3x10 ⁶	5.5	8/10

Washer/Disinfector

Removal/inactivation of inoculum (Non-Exposed) on Instruments

WD Conditions	Organism	Inoculum	Log Reduction	Positives
Routine	MRSA	2.6x10 ⁷	Complete	0/8
Routine	VRE	2.9x10 ⁷	Complete	0/8
Routine	P aeruginosa	2.1x10 ⁷	Complete	0/8
Routine	M terrae	1.2x10 ⁸	7.6	6/8
Routine	GS spores	8.1x10 ⁶	~1	12/12
No Enz/Det	VRE	2.4x10 ⁷	Complete	0/10
No Enz/Det	GS spores	8.7x10 ⁶	1.6	10/10

Washer/disinfectors are very effective in removing/inactivating microorganisms from instruments

Recommendations Monitoring of Sterilizers

- Monitor each load with physical and chemical (internal and external) indicators. If the internal indicator is visible, an external indicator is not needed.
- Use biological indicators to monitor effectiveness of sterilizers at least weekly with spores intended for the type of sterilizer (Class 6 emulating indicators not a substitute).
- Use biological indicators for every load containing implantable items and quarantine items, whenever possible, until the biological indicator is negative.

Types of Sterilization Monitoring Devices

- · Chemical Indicators
 - External chemical indicators
 - ◆Class 1 (process indicator, indicator tape)-outside of every package
 - Internal chemical indicators
 - ◆Class 2 (Bowie Dick)-routine testing of vacuum; within a test pack daily in an empty sterilizer
 - ◆Class 3 (single variable indicator; temperature, ETO conc)-may be used as internal monitor
 - ◆Class 4 (multi-variable indicator)-may be use as internal monitor

Types of Sterilization Monitoring Devices

- · Chemical Indicators
 - Internal chemical indicator
 - ◆Class 5 (integrating indicator)-may be used as internal monitor, suppose to mirnic the behavior of a biological indicator (BI)
 - ◆Class 6 (emulating indicator)-suppose to emulate or mirriic the behavior of a biological indicator, are cycle-specific (need a emulating indicator designed to valid a 10 min/270F cycle and a different indicator to validate a 3 min/270F). No professional organization (e.g., AORN, AAMI) has recommended the use of Class 6 emulating indicator as a substitute for biological indicators and there are no data that demonstrate that it mimics a B1 at suboptimal sterilization times.

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- · Current issues
 - Critical-cleaning with washer disinfectors, Class 6 chemical indicator
 - Semicritical items (e.g., endoscopes and C. difficile spores, new AERs)
 - Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores
 - → Microfiber
 - · Computers-sustained antimicrobial activity, touchscreen cleaning
 - ◆ Germicides-MRSA inactivation by disinfectants, technique
 - ♦ Green products

C. difficile spores

Disinfectants and Antiseptics C. difficile spores at 10 and 20 min, Rutala et al, 2006

- ~4 log₁₀ reduction (3 C. difficile strains including BI-9)
 - Clorox, 1:10, ~6,000 ppm chlorine (but not 1:50, ~1,200 ppm)
 - Clorox Clean-up, ~1,910 ppm chlorine
 - Tilex, ~25,000 ppm chlorine
 - Steris 20 sterilant, 0.35% peracetic acid
 - Cidex, 2.4% glutaraldehyde
 - Cidex-OPA, 0.55% OPA
 - Wavicide, 2.65% glutaraldehyde
 - Aldahol, 3.4% glutaraldehyde and 26% alcohol

High-Level Disinfection

C. difficile spores

- 2% glutaraldehyde is effective against C. difficile at 20
- 0.55% ortho-phthalaldehyde is effective against C. difficile at 10 minutes
- Steris 20 is effective against C. difficile at 10 and 20 minutes

Semicritical Equipment

- · Reprocessing semicritical items has been shown to have a narrow margin of safety
- Generally, the narrow margin of safety attributed to high microbial load and complex instruments with lumens
- · Any deviation from the recommended reprocessing protocol can lead to the survival of microorganisms and an increased risk of
- · Problems encountered with reprocessing semicritical equipment often related to improper cleaning

Errors in designing and reprocessing semicritical items continue and place patients at risk of infection

Automatic Endoscope Reprocessors (AERs)

- · Manual cleaning of endoscopes is prone to error.
- AER Advantages: automate and standardize reprocessing steps, reduce personnel exposure to chemicals, filtered tap water
- AER Disadvantages: failure of AERs linked to outbreaks, does not eliminate precleaning, does not monitor HLD concentration
- Problems: incompatible AER (side-viewing duodenoscope); biofilm buildup; contaminated AER; inadequate channel connectors; used wrong set-up or connector MMWR 1999;48:557
- · Must ensure exposure of internal surfaces with HLD/sterilant

Automatic Endoscope Reprocessors

- EvoTech-integrates cleaning (FDA-cleared claim) and disinfection.
 Automated cleaning comparable to manual cleaning. All residual data for cleaning of the internal channels as well as external insertion tube surfaces were below the limit of <8.5ug/cm²
- Reliance-requires a minimal number of connections to the endoscope channels and uses a control boot (housing apparatus the creates pressure differentials to ensure connectorless fluid flow through all channels that are accessible through the endoscope's control handle channel ports). Data demonstrate that the soil and microbial removal effected by Reliance washing phase was equivalent to that achieved by optimal manual cleaning. Alfa, Olson, DeGagne. AJIC 2006;34:561.

Infrared Coagulation (IRC)

- IRC is a widely used method for treating hemorrhoids. The procedure involves applying infrared light to compress and seal hemorrhoid veins.
- The manufacture sells a sterile disposable sheath and states removing and soaking lightguides between procedures is no longer required.
- The manufacture also states that the lightguide is damaged by immersion in a disinfectant (as the lightguide is not sealed at the end and disinfectant gets between the quartz glass and the covering)

Infrared Coagulation (IRC)

- CDC guideline (In press) recommends immersion for reprocessing endocavitary probes with covers because integrity of the cover is compromised (failure rate 1-80%)
- Since the lightguide cannot be immersed we investigated an alternative procedure
 - Wipe the probe for 2 minutes with 1:10 bleach
 - Wipe probe with sterile water and let air dry

Infrared Coagulation Testing

-:	· · · · · · · · · · · · · · · · · · ·	- Total Control of the Control of th		=
	Test Organism	Inoculum	Log ₁₀ Reduction (%)	
	Mycobacterium terrae	7.8 x 10 ⁶	6.9	

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- Current issues
 - Critical-cleaning with washer disinfectors, Class 6 chemical indicator
 - Semicritical items (e.g., endoscopes and *C. difficile* spores, new AERs)
 - Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores

 ◆ Wirrofiber
 - ◆ Computers-sustained antimicrobial activity, truchscreen cleaning
 - ◆ Germicides-MRSA inactivation by disinfectants, technique
 - ◆ Green products

Low-Level Disinfection for "Noncritical" Objects

Exposure time ≥ 1 min

Exposure time 2 1 min

Germicide Use Concentration

Ethyl or isopropyl alcohol 70-90%

Chlorine 100ppm (1:500 dilution)

Phenolic UD

lodophor UD

Quaternary ammonium UD

UD=Manufacturer's recommended use dilution

Disinfection and Sterilization of Emerging Pathogens

- · Hepatitis C virus
- Clostridium difficile
- Cryptosporidium
- Helicobacter pylori
- E.coli 0157:H7
- Antibiotic-resistant microbes (MDR-TB, VRE, MRSA)
- SARS Coronavirus, avian influenza, norovirus, prions
- Bioterrorism agents (anthrax, plague, smallpox)

C. difficile and Norovirus

Due to the relative resistance of *C. difficile* spores and norovirus, during clusters, surfaces should be disinfected with a product shown to be effective (e.g., chlorine 5000ppm [1:10 bleach]

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- · Current issues
- Critical-cleaning with washer disinfectors, Class 6 chemical indicator
- Semicritical items (e.g., endoscopes and C. difficile spores, new AERs)
- Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores
 - ◆ Microfiber
 - ◆ Computers-sustained antimicrobial activity, touchscreen cleaning
 - ◆ Germicides-MRSA inactivation by disinfectants, technique
 - ◆ Green products

Microfiber Cleaning

- Pad contains fibers (polyester and polyamide) that provide a cleaning surface 40 times greater than conventional string mops
- Proposed advantages: reduce chemical use and disposal (disinfectant solution not changed after every third room, clean microfiber per room [washing lifetime 500-1000x]); light (~5 lb less than string mop) and ergonomic; reduce cleaning times.
- Does the microfiber provide the same or better removal of microorganisms on surfaces?

Effectiveness of Microfiber Mop

- Test conditions with a EPA-registered disinfectant: compared routine mop and bucket; microfiber mop and bucket; microfiber mop and system bucket. Twenty-four replicates per condition.
- Conducted RODAC sampling before and after floor disinfection (5 samples per room)
- New disinfectant solution for each test condition
- Dry time varied from 2 (routine mop and bucket)-8 (microfiber mop and bucket) minutes

Effectiveness of Microfiber Mop (Rutale, Gergen and Weber, Am J Infect Control, 2007;35:569) Disinfectant-regular mop 95% Disinfectant-microfiber system 95% Disinfectant-microfiber mop and regular 88% mop bucket Detergent-regular mop 68%

95%

78%

Detergent-microfiber system

mop bucket

Detergent-microfiber mop and regular

Microfiber Summary

- The microfiber system demonstrated superior microbial removal compared to cotton string mops when used with a detergent cleaner
- The use of a disinfectant did not improve the microbial elimination demonstrated by the microfiber system
- Use of a disinfectant did significantly improve microbial removal when a cotton string mop was used

Disinfection of Computer Keyboards Computer Keyboards, ICHE 2006;27:372

- Increased use of computers in patient areas has led to contamination of keyboards as reservoirs of pathogens
- · Study performed to
 - Examine the efficacy of different disinfectants on the computer keyboard
 - Determine if there were cosmetic (key lettering removed) or functional changes after 300 wipes

Disinfection of Computer Keyboards

- All tested products were effective (>95%) in removing and/or inactivating the test pathogens (MRSA, P. aeruginosa). No functional/cosmetic damage after 300 wines.
- Disinfectants included: 3 quaternary ammonium compounds, 70% isopropyl atcohol, phenolic, chlorine (80ppm)
- At present, recommend that keyboards be disinfected daily (for 5 sec) and when visibly soiled

тавье 3.	Sustained I	Efficacy of	Disinfectants	Applied	o Keyboai	d Against	Vancomycin-Resistant
Enterococci	us Species					_	

	Effi		nfectant, by 'ation of Disi			nge
		nge at ours		nge at lours	Challe 48 I	nge at lours
Disinfectant	10-min Exposure	60-min Exposure	10-min Exposure	60-min Exposure	10-min Exposure	60-min Exposure
Alcohol	3.05	5.67	12.58	3.31	10.89	5.59
CaviWipes	100.00	100.00	100.00	100.00	100.00	100.00
Clorox Disinfecting Wipes	100.00	100.00	100.00-	100.00	100.00	100.00
Sani-Cloth Plus	100.00	100.00	100.00	100.00	100.00	100.00
Sterile water	0.00	0.28	9.69	0.00	0.00	9.09

NOTE. Efficiely was inheliated as the percentage difference in the number of colony-furming units on the treated keys, compared with the number of colony-forming units on the control keys, Challenge times are hours since disinfection exposure. QUATS demonstrated excellent sustained activity against VRE and antimicrobial activity was maintained over the 48 test period

Touchscreen Cleaning

- · Follow the manufacturer's recommendations
- · Prepare the cleaning solution according to the manufacturer's instructions (e.g., alcohol, glutaraldehyde, mild soap, phenolic)
- · Wet a clean, soft cloth with the selected cleaning solution
- · Remove excess liquid from the cloth and squeeze damp
- · Wipe exposed surfaces (do not allow liquid to enter interior)
- · Remove any soap residue by gently wiping with clean cloth
- · QUATS are not recommended by some manufacturers

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- Current issues
 - Critical-cleaning with washer disinfectors, Class 6 chemical indicator
 - Semicritical items (e.g., endoscopes and C. difficile spores, new AERs)
 - Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores
 - ◆ Microfiber
 - ◆ Computers-sustained antimicrobial activity, touchscreen cleaning
 - ◆ Germicides-MRSA inactivation by disinfectants, technique
 - ♦ Green products

MRSA

MRSA

- Frequency of environmental contamination in areas housing MRSA patients has ranged from 1 to 74% (23.1%, 53.6% from isolation rooms) of surfaces cultured.
- · MRSA viable in the environment for days to weeks
- . HCW can contaminate their hands or gloves by touching contaminated
- Cleaning or disinfecting the environment can reduce transmission but cleaning regimens, as currently practiced, may not eliminate MRSA from
- Since MRSA sensitive to all germicides, likely due to surfaces not cleaned/disinfected
- · Need targeted methods to evaluate the thoroughness of room cleaning

Practice or Product

Susceptibility of MSSA and MRSA to a Phenolic and Quaternary

Rutala et al. ICHE 1997;18:417

	Phenolic	Phenolic	QUAT	QUAT
	1:256	1:128	1:64	1:32
MSSA	2/60	0/60	5/60	1/60
MRSA	0/60	0/60	4/60	1/60

				enottous			
VS		VRd		MS:		MR	9A
0.5 ada	5 min	Q.5 min	6 astra	0.5 mln	5 rein	O.S mates	5 min
>4.3	×4.2	×4.8	>L3	>5.I	>5.I	>4.5	>4.0
NS.4	>5.4	***	>4.0	>5.0	>6.0	F-1.6	>4.5
×13	>4.3	>4.8	>4.S	95.5	>5.1	وترهد	>4.6
×5.5	>5.3	M\$.8	r6.5	×5.2	*b.1	×4.8	⊳4.6
9.1	5.3	1.9	3.7	+1.1	+0.9	-0.6	253
	9.5 ada 943 95.4 94.5 96.5	9.5 ards 5 min 943 942 954 954 954 955 943 943 943	0.6 ada 5 mia 0.5 mia 54.3 54.2 54.6 55.4 56.4 54.9 54.3 54.3 54.8 55.5 56.5 56.5 56.6	0.5 eds 5 mts 0.5 mts 6 mts >+3	0.5 mbs 5 mbs 0.5 mbs 0.5 mbs 1.5	0.5 mile 5 mile 0.3 mile 6 mile 0.5 mile T mile	0.5 ads 5 mds 0.5 mds 6 mds 0.5 mds 5 mds 0.5 mds 143 42 448 443 451 551 451 445 154 454 454 452 453 551 551 451 154 45 453 453 551 551 454 154 45 455 455 455 455 451 455 155 455 455 455 455 455 451 455

Not Product: Is It Practice?

Surface Disinfection Effectiveness of Different Methods

Technique (with cotton)	MRSA Log ₁₀ Reduction (QUAT)
Saturated cloth	4.56
Spray (10s) and wipe	4.56
Spray, wipe, spray (1m), wipe	4.56
Spray	4.56
Spray, wipe, spray (until dry)	4.56
Control: detergent	2.83

Patient Area Cleaning/Disinfecting PC Carling et al, SHEA 2007 and ICHE 2008;29:1

- · Monitor cleaning performance using an invisible fluorescent targeting method. Rooms (14 high-touch objects) were marked and evaluated after terminal cleaning.
- Results: 1,119 rooms and 13,369 objects were evaluated in 23 hospitals. Mean proportion of objects cleaned was 49%. Following education and process improvement feedback, cleaning improved
- · Conclusion: Substantial opportunity for improving terminal cleaning/disinfecting activities.

Practice* NOT Product

*surfaces not wiped

The Green Hospital Keslar. Proto (Mass Gen Hosp). Fall 2007

- Hospitals are feeling the pressure to go green, both from eco-conscious donors and governmental agencies
- · Some features of The Green Hospital
 - Roof garden-wildlife habitat, adds insulation, absorbs rain
 - Fewer contaminants-uphoistery and mattresses without flame retardants, formaldehyde-free insulation, green cleaning products (no hazardous furnes), triple-filtered air
 - Exposure to natural light

 - Reduced water usage-water efficient toilets and faucets
 Greater energy efficiency-low energy fluorescent bulbs
 - More quiet-number 1 complaint is noise, better insulation between rooms

				Log . Re	ductions			
	Stapbylo		Salmo choicta	sella	Escher co# 015		Psetido. aerugi	
Product	0,5 min	5 min	0.5 mln	S walk	9.5 ann	5 เพล	0.5 min	5 m m
Vesphene II.se	58.2	×8.2	>6.7	×6.7	પ્રદેશ	>5.6	>6.7	3 N T
TEQ	+6.4	46.4	≥6.5	v6:5	>6.4	>6.4	×6.9	>62
Closox	55.S	ან.გ	5S.9	25.9	>\$.6	:5.6	2 5 .3	55.2
Ethanel	5.2	>6.7	>6.0	×6.0	>6.6	>6.8	>6.1	>6.4
Lysol Disinfectant	12	4.3	4.0	19	40	4.1	4.2	4.0
Lysol Antikecterial	>5.6	+56	35.8	>5.5	>5.7	>5.7	55.5	×5.5
Mr. Clean	4.1	+6.0	>5.7	-5.7	×6.1	4.7	>5.7	-5.7
Vinegar	9.03	6.3	56,0	56B	0.4	2.4	s5.6	35.6
Baking seda	0.2	0.5	2.3	2.3	6.4	6.7	1.6	1.1

			X 110CL 7C-0 1C-		STANT BACTER			
	VSE		VRE		MSSA		MRSA	
Product	0.5 mta	C min	Q.6 min	5 onta	0.5 mla	15 mm	0,5 min	F noise
Vesobene ilse	~4.3	>4.2	અહ	24.3	>5.1	25.1	×4.6	24.6
Clorux	>5-4	×5.4	>4.0	143	>6.0	>5.0	>4.6	>4.3
Lysel Disinfectant	×3	>4.2	4.3	244	25.5c	-S-1	>4.6	÷.5
Lysot Autibucterial	>5.5	25.5	5-5.5	>53	≥5.1	1.6<	>4.6	>4.4
Vinerar	0.1	5.3	1.0	3.7	+1.1	-0.0	+0.6	2.3

Rutata WA, Barbee SL, Aguiar NC, Sobsey MD, Weber DJ. Antimicrobial Activity of Home Disinfectants and Natural Products Against Potential Human Pathogens, Infection Control and Hospital Epidemiology 2000;21:33-38

"Green" cleaners will remove microbial contaminants but will not dependably kill

Summary

- D/S guidelines must be followed to prevent exposure to pathogens that may lead to infection. Semicritical items represent the greatest risk. Class 6 indicators not a substitute for biological indicators.
- During clusters, surfaces potentially contaminated with norovirus or C. difficile spores should be disinfected with with an agent shown to have efficacy (e.g., hypochlorite, 5000 ppm)
- Washer-disinfectors are very effective in removing/inactivating microbes
- Microfiber demonstrated superior microbial removal compared to cotton-string mops with a detergent
- Disinfectants (but not natural products) demonstrate excellent activity against MRSA but practices are deficient. QUATS have sustained antimicrobial activity.

Current Issues in Disinfection and Sterilization

- · Disinfection and sterilization principles
- Current issues
 - Critical-cleaning with washer disinfectors, Class 6 chemical indicator
 - Semicritical items (e.g., endoscopes and C. difficile spores, AERs)
 - Noncritical-surface disinfection
 - ◆ Norovirus and C. difficile spores
 - ◆ Microfiber
 - ◆ Computers-sustained antimicrobial activity, touchscreen cleaning
 - ◆ Germicides-MRSA inactivation by disinfectants, technique
 - ◆ Green products

disinfectionandsterilization.org

Thank you

References

- Rutala WA, Weber DJ. CJD: Recommendations for disinfection and sterilization. Clin Infect Dis 2001;32:1348
- Rutala WA, Gergen MF, Weber DJ. Microbiologic evaluation of microfiber mops for surface disinfection. Am J Infect Control 2007;35:569
- Rutala WA, Weber DJ, HICPAC. CDC guideline for disinfection and sterilization in healthcare facilities. MMWR. In press.
- Rutala WA. APIC guideline for selection and use of disinfectants. Am J Infect Control 1996;24:313

References

- Rutala WA, Peacock JE, Gergen MF, Sobsey MD, Weber DJ. Efficacy of hospital germicides against adenovirus 8, a common cause of epidemic keraloconjunctivitis in health care facilities. Antimicrob Agents Chemother 2006;50:1419
- Rutala WA, White MS, Gergen MF, Weber DJ. Bacterial contamination of keyboards: Efficacy and functional impact of disinfectants. Infect Control Hosp Epidemiol 2006;27:372
- Rutala WA, Weber DJ. Surface disinfection: Should we do it? J Hosp Infect. 2000; 48:S64.
- Schneider PM. New technologies for disinfection and sterilization. In: Rutala WA (ed). Disinfection, Sterilization and Antisepsis. 2004:127-139